

Fluid processor

Booth 2511 - Designed to operate at up to 30,000 psi process pressure, the M110EH-30 laboratory Microfluidizer processor, is beneficial to the production of nano-suspensions and nano-emulsions as well as cell disruption with the least number of passes. Additionally, more efficient processing of pharmaceutical formulations, fine chemicals, and biological material ranging from simple oil-in-water emulsions to highly immiscible solids-in-liquid suspensions will result.

Source URL (retrieved on 01/29/2015 - 3:14am):

http://www.chem.info/news/2006/03/fluid-processor?qt-recent_content=0&qt-most_popular=1