

Data Management for pH & Ion Measurements



Mettler Toledo has released its LabX® Direct pH software for data management of pH and ion measurements. The software features:

- No manual transcription, paper printouts or copying and pasting data manually from one application to another.
- Flexible data transfer, which sends data from any Mettler Toledo SevenGo™, SevenEasy™ or SevenMulti™ meter (Version 2.0) to Microsoft® Excel® or any other Windows-based application.
- Flexible transfer parameters, which can be selected for transfer to an Excel spreadsheet.
- Support for the connection of a barcode reader to help reduce errors and speed up tasks.
- Capability to display measurement curves.

www.mt.com [1]

Source URL (retrieved on 03/06/2015 - 4:47am):

Data Management for pH & Ion Measurements

Published on Chem.Info (<http://www.chem.info>)

<http://www.chem.info/product-releases/2010/01/data-management-ph-ion-measurements>

Links:

[1] <http://www.mt.com>