

Filter Facilitates Single-Use Integration



Meissner has designed its CM and CK capsule filters for ease of integration into single-use systems. According to the company, the capsule filter further delivers:

- The capability to minimize fluid hold-up volume for maximum product recovery.
- Compatibility with small-volume, single-use biocontainers and tubing assemblies.
- A full range of Meissner membranes, including PES, PVDF, PTFE and PP.
- Availability of no vents, 1 vent or 2 vents in either luer ports or sanitary vent valves.
- The flexibility to be specified with a range of PP microfiber media and glass fiber.
- A trim profile.
- Optional gamma-irradiation.

leads@meissner.com [1]

www.meissner.com [2]

Source URL (retrieved on 01/29/2015 - 8:51am):

<http://www.chem.info/product-releases/2012/04/filter-facilitates-single-use->

Filter Facilitates Single-Use Integration

Published on Chem.Info (<http://www.chem.info>)

[integration?cmpid=related_content](#)

Links:

[1] <mailto:leads@meissner.com>

[2] <http://www.meissner.com/>