

## Do without Troublesome Liners



VICI Valco has designed its 360-micron high-pressure fittings to permit the direct connection of OD fused silica, PEEK, stainless or electroformed nickel tubing without having to use troublesome liners. According to the company, other features include:

- A leak-free connection, which can be easily be formed with the available manual tool, that seals at pressures well in excess of 20,000 PSI.
- A ferrule that is free to rotate and does not twist the tube as the nut is tightened, as well as that snaps into the nut, so that the fitting is 1 piece.
- A compact size and fine 2-56 threads.
- A wide range of micron bore sizes.
- Availability with a PEEK or stainless nut; the PEEK nut is supplied with a glass-filled PEEK ferrule, while the stainless nut uses a 316 stainless steel ferrule.
- Suitability for laboratories, food and beverage, pharmaceutical and chemical applications.

[valco@vici.com](mailto:valco@vici.com) [1]

[www.vici.com](http://www.vici.com) [2]

**Source URL (retrieved on 02/01/2015 - 1:40pm):**

---

## **Do without Troublesome Liners**

Published on Chem.Info (<http://www.chem.info>)

---

[http://www.chem.info/product-releases/2012/06/do-without-troublesome-liners?qt-most\\_popular=1](http://www.chem.info/product-releases/2012/06/do-without-troublesome-liners?qt-most_popular=1)

### **Links:**

[1] <mailto:valco@vici.com>

[2] <http://www.vici.com/>