

## Taking More than One Position



Torrey Pines Scientific presents its 2-position *EchoTherm*™ IC22XT digital chilling and heating dry bath for freezing, chilling or heating biological samples from -20° to 100°C. According to the company, the *EchoTherm* is additionally equipped with:

- The ability to run 2 separate temperatures and 2 different sample blocks at the same time.
- A built-in data logger, as well as a 30-day countdown timer with an alarm.
- An RS232 interface, making it suitable for use with robotic systems.
- 2 dry baths in 1 to save bench space.
- Models with 100, 115 or 230 VAC, 50/60 Hz.
- Suitability for ligations at 14°C; maintaining 17°C for storing oocytes; enzyme reactions and deactivations; hybridizations; incubating samples above, below or at room temperature; storing enzymes or DNA libraries, storing samples at ice bucket temperatures; and more.
- Compliance with UL, CSA and CE or equivalent agency certifications.
- The flexibility to be used with accessory sample blocks for 0.5-, 1.5-, and 2.0-mL centrifuge tubes, PCR tubes and plates, 96-well and 384-well assay plates of all kinds, vials and most sizes of test tubes.
- Availability of custom blocks.

## Taking More than One Position

Published on Chem.Info (<http://www.chem.info>)

---

[www.torreypinesscientific.com](http://www.torreypinesscientific.com) [2]

### Source URL (retrieved on 02/01/2015 - 2:15pm):

[http://www.chem.info/product-releases/2012/10/taking-more-one-position?qt-most\\_popular=0](http://www.chem.info/product-releases/2012/10/taking-more-one-position?qt-most_popular=0)

### Links:

[1] <mailto:info@torreypinesscientific.com>

[2] [http://app.bronto.com/public/?q=ulink&fn=Link&ssid=7172&id=if9j7x9irj90hppjxpvz6ubiabmhh&id2=c16c5kr4ipx2p5ybfre5n8r6rldpe&subscriber\\_id=bofzzzdmydbxcrsabstuxxcrylrbnm&delivery\\_id=aulryzasjgngbrkkbzbxbatsdbtnban&tid=3.HAQ.ALL0mQ.FR08.QcLm..djlH.b..l.BOH9.a.UCKpIA.UCKzrA.Oskwlw](http://app.bronto.com/public/?q=ulink&fn=Link&ssid=7172&id=if9j7x9irj90hppjxpvz6ubiabmhh&id2=c16c5kr4ipx2p5ybfre5n8r6rldpe&subscriber_id=bofzzzdmydbxcrsabstuxxcrylrbnm&delivery_id=aulryzasjgngbrkkbzbxbatsdbtnban&tid=3.HAQ.ALL0mQ.FR08.QcLm..djlH.b..l.BOH9.a.UCKpIA.UCKzrA.Oskwlw)